


Clinical Medical Image Library

Image 1:007

Drusen on Fire: Ex Vivo Hyperspectral Autofluorescence (AF) Image of Drusen from 81-Year-Old Male Donor with Age-Related Macular Degeneration



The image is an overlay of 3 images representing the distribution of the 3 most abundant fluorophore signals that were extracted using advanced factorization algorithms. The 3 images were digitally colored red (1), green (2) and blue (3). The signals localized to different components in the RPE as illustrated by color mixing: normal RPE cells (light blue), RPE Cells stacked at the edge of drusen creating a halo effect (light green), main body of drusen (blue) with a distinct central signal (pink) that also localizes to sub-RPE deposits (red). This method demonstrates that distinct fluorophore signals are localized specifically to drusen in AMD tissues.

Information

R. Theodore Smith*, Tal Ben-Ami and Yuehong Tong

New York University, Department of Ophthalmology, Newyork, USA

*Correspondence: R. Theodore Smith, New York University, Department of Ophthalmology, 550 First Avenue, New York, NY 10016, USA, Tel: 212-263-7300; E-mail: tbenami@gmail.com

Citation: Smith RT, Ben-Ami T, Tong Y (2015) Drusen on Fire: Ex Vivo Hyperspectral Autofluorescence (AF) Image of Drusen from 81-Year-Old Male Donor with Age-Related Macular Degeneration. Clin Med Img Lib 1:007

Published: August 25, 2015

Copyright: © 2015 Smith RT. This is an open-access content distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

